

Sômatik Labs®



SLS-IRCT30

Handheld Infrared
Contactless Thermometer

USER MANUAL

CONTENTS

1.0	BEFORE USE	3
2.0	SAFETY PRECAUTIONS	3
3.0	PACKAGE CONTENTS	4
4.0	NORMAL BODY READING TEMPERATURE.....	4
4.1	Variation in body temperature.....	5
5.0	PRODUCT DESCRIPTION	5
5.1	Overview.....	5
5.2	Operating principle	5
5.3	Intended use.....	5
5.4	Features.....	6
6.0	PRODUCT OVERVIEW.....	7
6.1	Interface overview	7
7.0	SOUND & LCD COLOUR INDICATOR	8
7.1	Forehead temperature.....	8
7.2	Object temperature	8
8.0	DISPLAY AND OPERATING INSTRUCTIONS	9
8.1	Measuring Body temperature.....	9
8.2	Measuring Object temperature	9
8.3	Out of the measuring range display	10
8.4	Switching between °C and °F.....	10
8.5	Switching between forehead temperature & object temperature.....	11
8.6	Switching between mute and un-mute.....	11
8.7	Recall memory	11
8.8	Change default fever alert threshold settings ..	12
8.9	Error information & low battery	13
8.10	Measurement Process.....	13

9.0	REPLACING BATTERIES	14
10.0	CLEANING	15
11.0	MAINTENANCE	15
12.0	TROUBLESHOOTING	16
13.0	SPECIFICATIONS.....	17
13.1	Guidance and manufacturer's declaration - Electromagnetic emission - for all equipment and systems	18
13.2	Guidance and manufacturer's declaration - Electromagnetic immunity - for all equipment and systems	19
13.3	Guidance and manufacturer's declaration - Electromagnetic immunity - for equipment and systems that are not life supporting	20
13.4	Recommended separation distances between portable and mobile RF communications equipment and the equipment or system - for equipment and systems that are not life supporting	22
14.0	WARRANTY TERMS & CONDITIONS	23
14.1	Indemnity.....	24
14.2	Warranty and Warranty Disclaimer	24
15.0	NOTES	28

1.0 BEFORE USE

The Somatik Labs SLS-IRCT30 Contactless Thermometer has been carefully developed for accurate, safe and fast temperature measurements. As with any thermometer it is important to follow the included instructions to ensure accurate temperature readings. Please open and check the contents carefully before use and perform all installation and operating instructions as per this user manual. In case of any damage or operation problem please contact Somatik Labs service department (refer page 23).

2.0 SAFETY PRECAUTIONS

Read the following precautions carefully before using the thermometer.

ATTENTION

- Thermometer Probe is not be touched or tampered with.
- To protect the environment, please dispose of the batteries responsibly. Warning: Batteries may explode is disposed of in fire.
- Please remove the batteries if you don't intend to use the thermometer for more than 2 months.
- Do not immerse the thermometer in water.
- Do not subject the thermometer to vibration, shock or impact.
- Do not take body temperature readings within 20 minutes after you do physical exercises as the reading may not be accurate of your body temperature.
- Do not touch the tip of the temperature probe as it is directly connected with the sensor.
- In case the ambient temperature is either too hot or too cold, please let the thermometer adjust to room temperature by keeping it idle for 30 minutes before use.
- Operating temperature of the thermometer is between 10°C (50°F) and 40°C (104°F).
- Thermometer is only compatible with 1.5V AAA batteries. Please do not use the batteries with different voltages or specifications.

WARNINGS & DISCLAIMERS

- The thermometer is not intended to diagnose or treat any health problem, measurement results are for reference only.
- It is dangerous to make a self-diagnosis or self-treatment based on the obtained measurement results. For such purposes, please consult a physician or other medical professionals.
- Do not disassemble the thermometer or attempt to repair it, thermometer may get damaged permanently.
- During measurement, do not use a mobile phone or any other device that may cause electromagnetic interference.
- Please keep the thermometer out of the reach of children.

3.0 PACKAGE CONTENTS

1. Infrared Thermometer x 1
2. Soft Pouch x 1
3. Battery (AAA) x 2
4. Instruction Manual x 1

4.0 NORMAL BODY READING TEMPERATURE

Body Part:

Forehead


Normal Temperature Range:

36.1°C–37.5°C / 97.0°F–99.5°F

The normal body temperature range slightly varies with age and gender. Generally, newborns or children have higher body temperature than adults.

Women's temperature is approximately 0.3°C higher than men's.

4.1 VARIATION IN BODY TEMPERATURE



Normal body temperature varies by the time of day and is also affected by external factors.

The body temperature of an individual is the lowest between 2:00 a.m. and 4:00 a.m. and the highest between 14:00 p.m. and 20:00 p.m. An individual's body temperature typically changes by less than 1°C each day.

5.0 PRODUCT DESCRIPTION

Infrared Thermometer SLS-IRCT30 measures the human body's or an object's temperature based on the infrared energy emitted by the forehead or an object (such as milk and water). You can quickly get measurement results after pointing the temperature probe at the target and then press the trigger.

5.1 OVERVIEW

The thermometer consists of plastic body, LCD to display, buttons, buzzer, infrared temperature sensor and microprocessor circuit.

5.2 OPERATING PRINCIPLE

The infrared temperature sensor collects infrared energy emitted by the forehead. The energy is converted into a temperature reading by the thermocouple and the measurement circuit.

5.3 INTENDED USE

The Infrared Thermometer SLS-IRCT30 is a non-contact infrared thermometer intended to obtain the body temperature from the forehead or of object.

5.4 FEATURES

SAFE TO OPERATE

- Passive infrared receiving technology.
- Non-contact measurement, to prevent cross-infection.

EASY OPERATION

- Handheld design, easy operation.
- One-click automatic temperature measurement.

QUICK RESPONSE

- 1 second measurement.


HIGH ACCURACY

- Advanced infrared temperature sensor, with high sensitivity.
- Enhanced accuracy with automatic temperature calibration.

FUNCTIONS

- 20 temperature readings stored in memory.
- Forehead/Object temperature measurement.
- Fever alert, with a configurable alert threshold.
- Switching between °C and °F.
- Mute/Un-Mute mode (measuring sound notification).
- Auto power-off.


6.0 PRODUCT OVERVIEW



1. IR sensor
2. Power-on trigger / Measure trigger
3. Battery cover

4. LCD display
5. Mute / Un-mute button
6. Mode button
7. Celsius / Fahrenheit Switch button
8. Handgrip

6.1 INTERFACE OVERVIEW



1. Object temperature mode
2. Body temperature mode
3. Temperature unit ($^{\circ}\text{C}$ / $^{\circ}\text{F}$)
4. Low power indicator
5. Mute / un-mute
6. Memory mode
7. Temperature reading

7.0 SOUND & LCD COLOUR INDICATOR

7.1 FOREHEAD TEMPERATURE

RANGE	SOUND	SCREEN COLOUR
34.9°C - 37.5°C 94.8°F - 99.5°F	1 long beep	Green
37.6°C - 42.2°C 99.6°F - 108.0°F	3 short double beeps	Red

7.2 OBJECT TEMPERATURE


RANGE	SOUND	SCREEN COLOUR
0°C - 100°C 32.0°F - 212.0°F	1 long beep	White

NOTE:


- When the temperature reading is between 34.9°C/94.8°F and 37.5°C/99.5°F, there will be a long beep and a green backlight which indicates normal temperature reading.
- When the temperature reading is between 37.6°C/99.6°F and 42.2°C/108.0°F, there will be 3 short double beeps and a red backlight. This indicates that the body temperature is a little high and you may have a fever. Please consult your doctor if you are not sure.

8.0 DISPLAY AND OPERATING INSTRUCTIONS



8.1 MEASURING BODY TEMPERATURE

Screen Display	Description	Audio & Visual Alert
	<p>In a power-off state, point the IR sensor to the center of the forehead. Move the thermometer towards the forehead.</p> <p>For effective measurement, the distance between the thermometer and the forehead must be 1-5 cm (½" to 2").</p> <p>Press and release the Trigger button. The forehead temperature will be displayed on the screen.</p>	See the table in the "Sound and LCD Colour Indicator" section



8.2 MEASURING OBJECT TEMPERATURE

Screen Display	Description	Audio & Visual Alert
	<p>In a power-on state, press the "Mode button", the thermometer enters the Object mode.</p> <p>Point the IR sensor to the center of the object, then press and release the Trigger button. The object temperature will be displayed on the screen.</p>	See the table in the "Sounds and Backlight Instructions" section





8.3 OUT OF THE MEASURING RANGE DISPLAY

Screen Display	Description	Audio & Visual Alert
	In Object mode, a temperature reading of more than 100°C (212.0°F)	A long beep with green LCD colour for 3 seconds.
	In Body Temperature mode, a temperature reading of more than 42.2°C (108.0°F)	
	In Object mode, a temperature reading of less than 0°C (32.0°F)	A long beep and a green LCD colour for 3 seconds.
	In body temperature mode, a temperature reading of less than 34.9°C (94.8°F)	





8.4 SWITCHING BETWEEN °C AND °F

Screen Display	Description	Audio & Visual Alert
 	In a power-on state, press the °C/°F button to switch between °C and °F.	Silent


8.5 SWITCHING BETWEEN FOREHEAD TEMPERATURE & OBJECT TEMPERATURE

Screen Display	Description	Audio & Visual Alert
 	In a power-on state, press the Mode button to switch between forehead temperature  and object temperature  .	Silent

8.6 SWITCHING BETWEEN MUTE AND UN-MUTE

Screen Display	Description	Audio & Visual Alert
 	In a power-on state, press the Mute/Un-mute button  to switch between mute and un-mute.	The  symbol is displayed in Un-mute mode and disappears in mute mode.

8.7 RECALL MEMORY




Screen Display	Description	Audio & Visual Alert
	In power-on state, press and hold the Mode button for more than 2 seconds. "F-1" is displayed.	Press the Measure button to return to the measurement interface.

Screen Display	Description	Audio & Visual Alert
	<p>Press the °C/°F or the button, 1 will be shown, followed by the recorded reading.</p>	Silent
	<p>Press the "°C/°F button" again for the next recorded data. 2 will be shown, followed by the recorded reading.</p>	
	<p>A maximum of 20 temperature readings can be recalled.</p>	
	<p>Note: 1 represents the newest data.</p>	

8.8 CHANGE DEFAULT FEVER ALERT THRESHOLD SETTINGS

Screen Display	Description	Audio & Visual Alert
	<p>When "F-1" is displayed, press the Mode button. Then "F-2" is displayed.</p>	<p>Press the Measure button to return to the measurement interface.</p>
	<p>Press the °C/°F or the button. The fever alert threshold is displayed. The threshold value increments by 0.1°C/°F every time the °C/°F button is pressed, and decrements by 0.1°C/°F every time the button is pressed. The setting range is 35.0°C–42.0°C (95.0°F–107.6°F).</p>	<p>The default fever alert threshold is changed to $\geq 37.6^{\circ}\text{C}$.</p>

8.9 ERROR INFORMATION & LOW BATTERY

Screen Display	Description	Audio & Visual Alert
	The ambient temperature is higher than 40.0°C (104.0°F) or lower than 10.0°C (50.0°F).	A long beep and a red backlight for 3 seconds.
	The recording was not complete.	A long beep and a red backlight for 3 seconds.
	Battery has low power, please replace the batteries.	Silent
Power-off	In any mode, if there is no operation for 10 seconds, the thermometer will power off automatically.	

8.10 MEASUREMENT PROCESS

1. Select the measurement mode.
 - Press the Measure button to power on the thermometer. Select the measurement mode using the Mode button.
 - The Ω symbol indicates the Body Temperature mode. The \square symbol indicates the Object Temperature mode.
2. Press the Measure button to start a measurement.
 - When taking the forehead temperature, point the IR sensor to the center of the forehead. Move the thermometer towards the forehead. The distance between the thermometer and the forehead must be 1-5 cm (½" to 2"). Press and release the Trigger button. The forehead temperature will be displayed on the screen.

- When taking the object temperature, Point the IR sensor to the center of the object. The distance between the thermometer and the object must be 1-5 cm (½" to 2"). Press and release the Trigger button. The object temperature will be displayed on the screen.



- After each measurement, please clean the thermometer with a dry soft cloth, and put the thermometer in a dry and well-ventilated place.
- The thermometer automatically powers off if it is not used in 10 seconds.

NOTES:

- The thermometer is suitable for an indoor environment without strong air convection between the thermometer and the target. For example, winds from a fan, an air-conditioner, or a heater.
- Do not hold the thermometer for a long time, because it is sensitive to the ambient temperature.
- Make sure the sensor head is free of foreign matters before use.
- Make sure the forehead has no sweat and no hairs covered before measure the forehead temperature; otherwise, the result could be incorrect.

9.0 REPLACING BATTERIES

1. At the base of the thermometer slide the battery cover off along the marked direction. Insert the two AAA batteries into the compartment correctly.
2. If the low-battery symbol is displayed on the screen, replace the batteries.

NOTES:

- Make sure that the batteries are installed correctly. Otherwise, the thermometer may be damaged.
- Batteries of a same type should be used. Dispose the used batteries in accordance with the local environmental policies.
- The thermometer is provided with batteries already installed. When you use it in the first time, please open the battery cover, and remove the insulation cover protecting the battery discharge.

10.0 CLEANING**CLEANING STEPS:**

1. Take the batteries out before cleaning.
2. Clean the temperature with a soft cloth. Clean the lens of the temperature probe with a cotton swab.
3. Wipe the thermometer body with a slightly damp soft cloth.

NOTES:

- The lens may be scratched if it is cleaned with a piece of tissue paper, which might result in inaccurate readings.
- Do not clean the thermometer with corrosive cleansers. During the cleaning process, do not touch the lens using hard objects, do not immerse any part of the thermometer into liquid, or allow liquid to enter into the thermometer.

11.0 MAINTENANCE

1. Ensure the safety of the thermometer and check prior to use for any damage to the lens or cracks in the housing. Do not use the thermometer if there is any obvious damage. Please store the thermometer in the provided storage bag when not in use.

2. Store the thermometer in a dry, dust-free, and well-ventilated place. Make sure that the thermometer is not exposed to sunlight. Make sure that the storage and transportation environments meet the requirements.
3. Remove the batteries if the thermometer will not be used for more than two months.

12.0 TROUBLESHOOTING

THERMOMETER FAILS TO POWER ON

Possible Cause:

1. Low battery.
2. Polarities of the batteries are reversed.

Solution:

1. Change the batteries.
2. Make sure that the batteries are installed correctly.

"ER1" IS DISPLAYED.

Possible Cause:

1. The ambient temperature is lower than 10°C (50.0°F) or higher than 40°C (104°F). Polarities of the batteries are reversed.

Solution:

1. Take a measurement under an ambient temperature between 10°C (50.0°F) and 40°C (104°F).

TEMPERATURE READING IS LOWER THAN THE TYPICAL BODY TEMPERATURE RANGE

Possible Cause:

1. The lens of the temperature probe is dirty.
2. The distance between the temperature probe and the target is too long.
3. The thermometer is used within 30 minutes after being taken from a cold environment.

Solution:

1. Clean the lens using a cotton swab.
2. Move the thermometer closer to the target.
3. Wait for more than 30 minutes after the thermometer is moved into the measurement environment.

TEMPERATURE READING IS HIGHER THAN THE TYPICAL BODY TEMPERATURE RANGE

Possible Cause:

1. The temperature probe is damaged.

Solution:

1. Please contact our service team.

13.0 SPECIFICATIONS

Product Name	Infrared Thermometer
Product Model	SLS-IRCT30
Operating Voltage	DC 3V
Battery Model	AAA x 2
Display	Segment LCD
Measure time	About 1 second
Measuring Distance	1 to 5cm (1/2" to 2")
Measuring Range	Forehead: 34.9°C–42.2°C (94.8°F–108.0°F) Object: 0.0°C–100.0°C (32.0°F–212.0°F)
Accuracy (Laboratory)	$\pm 0.4^{\circ}\text{F}/\pm 0.2^{\circ}\text{C}$ from 94.8°F to 108.0°F (34.9°C to 42.2°C) $\pm 0.5^{\circ}\text{F}/\pm 0.3^{\circ}\text{C}$, Outside the range of 94.8°F to 108.0°F (34.9°C to 42.2°C)
Resolution	0.1°C (0.1°F)
Memory	20 temperature readings
Low-battery Alert	The low-battery symbol is displayed if the power voltage is lower than 2.5 V \pm 0.1V.
Automatic Power-off	The thermometer automatically powers off if it is not used in 10 \pm 1 seconds.
Dimensions (mm)	150 \times 88.2 \times 40.6
Weight (g)	109.5 g (with batteries)

Operating Environment	Temperature: 10°C–40°C (50°F–104°F)
	Humidity: 15%–95% RH, non-condensing
	Atmospheric pressure: 86–106 kPa
Storage and Transportation	Temperature: -20°C to 55°C (-4°F–131°F)
	Humidity: 0–95% RH, non-condensing
	Atmospheric pressure: 50–106 kPa
<p>The infrared thermometer has been tested and conforms to the standard ASTM E1965-98. ASTM laboratory accuracy requirements in the display range of 98°F to 102°F (37°C–39°C) for skin IR thermometers is $\pm 0.5^\circ\text{F}$ ($\pm 0.3^\circ\text{C}$). Note that for mercury-in-glass and electronic thermometers, the requirement per ASTM Standards E667-86 and E1112-86 is $\pm 0.2^\circ\text{F}$ ($\pm 0.1^\circ\text{C}$).</p>	


13.1 GUIDANCE AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC EMISSION - FOR ALL EQUIPMENT AND SYSTEMS

Guidance and manufacturer's declaration - Electromagnetic immunity		
<p>The Infrared Thermometer SLS-IRCT30 is intended for use in the electromagnetic environment specified below. The customer or the user of the Infrared Thermometer SLS-IRCT30 should assure that it is used in such an environment.</p>		
Emissions test	RF emissions CISPR 11	RF emissions CISPR 11
Compliance	Group 1	Class B
Electromagnetic environment - guidance	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

13.2 GUIDANCE AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC IMMUNITY - FOR ALL EQUIPMENT AND SYSTEMS

Guidance and manufacturer's declaration - Electromagnetic immunity		
The Infrared Thermometer SLS-IRCT30 is intended for use in the electromagnetic environment specified below. The customer or the user of the Infrared Thermometer SLS-IRCT30 should assure that it is used in such an environment.		
Immunity test	Electrostatic discharge (ESD) IEC 61000-4-2	Power frequency (50/60 Hz) magnetic field IEC 61000-4-8
IEC 60601 test level	±6 kV contact ±8 kV air	3 A/m
Compliance level	±6 kV contact ±8 kV air	3 A/m
Electromagnetic environment - guidance	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

13.3 GUIDANCE AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC IMMUNITY - FOR EQUIPMENT AND SYSTEMS THAT ARE NOT LIFE SUPPORTING

Guidance and manufacturer's declaration - Electromagnetic immunity	
The Infrared Thermometer SLS-IRCT30 is intended for use in the electromagnetic environment specified below. The customer or the user of the Infrared Thermometer SLS-IRCT30 should assure that it is used in such an environment.	
Immunity test	Radiated RF IEC 61000-4-3
IEC 60601 test level	3 V/m 80 MHz to 2.5 GHz
Compliance level	3 V/m
Electromagnetic environment - guidance	<p>Portable and mobile RF communications equipment should be used no closer to any part of the SLS-IRCT30, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = \left[\frac{3.5}{E_1} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = \left[\frac{7}{E_1} \right] \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ <p>Where p is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).^b</p> <p>Field strengths from fixed RF transmitters, as determined by an Electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol: </p>

NOTE:

- At 80 MHz and 800 MHz, the higher frequency range applies.
- These guidelines may not apply in all situations. Electromagnetic is affected by absorption and reflection from structures, objects and people.

a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy.

To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the SLS-IRCT30 is used exceeds the applicable RF compliance level above, the SLS-IRCT30 should be observed to verify normal operation.

If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the SLS-IRCT30.

b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

13.4 RECOMMENDED SEPARATION DISTANCES BETWEEN PORTABLE AND MOBILE RF COMMUNICATIONS EQUIPMENT AND THE EQUIPMENT OR SYSTEM - FOR EQUIPMENT AND SYSTEMS THAT ARE NOT LIFE SUPPORTING

The Infrared Thermometer SLS-IRCT30 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Infrared Thermometer SLS-IRCT30 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Infrared Thermometer SLS-IRCT30 as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W.	Separation distance according to frequency of transmitter.	
	80 MHz to 800 MHz $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$
0.01	0.12	0.23
0.1	0.38	0.73
1	1.2	2.3
10	3.8	7.3
100	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE:

- At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.
- These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

14.0 WARRANTY TERMS & CONDITIONS

Our goods come with guarantees that cannot be excluded under the Australian Consumer Law. You are entitled to a replacement or refund for a major failure and for compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the goods repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure.

This warranty is provided in addition to your rights under the Australian Consumer Law. Directed Electronics warrants that this product is free from defects in material and workmanship for a period of 12 months from the date of purchase or for the period stated on the packaging. This warranty is only valid where you have used the product in accordance with any recommendations or instructions provided by Directed Electronics.

This warranty excludes defects resulting from alterations of the product, accident, misuse, abuse or neglect. In order to claim the warranty, you must return the product to the retailer from which it was purchased or if that retailer is part of a National network, a store within that chain, along with satisfactory proof of purchase. The retailer will then return the goods to Directed Electronics.

Directed Electronics will repair, replace or refurbish the product at its discretion. The retailer will contact you when the product is ready for collection. All costs involved in claiming this warranty, including the cost of the retailer sending the product to Directed Electronics, will be borne by you.

Email: service@somatiklabs.com

Phone: +61 03 8331 4800

14.1 INDEMNITY

YOU AGREE TO DEFEND, INDEMNIFY AND HOLD HARMLESS SOMATIK LABS AND ITS SUBSIDIARIES AND AFFILIATES FROM AND AGAINST ANY AND ALL CLAIMS, PROCEEDINGS, INJURIES, LIABILITIES, LOSSES, COSTS AND EXPENSES (INCLUDING REASONABLE LEGAL FEES), INCLUDING BUT NOT LIMITED TO, CLAIMS ALLEGING NEGLIGENCE, INVASION OF PRIVACY, COPYRIGHT INFRINGEMENT AND/OR TRADEMARK INFRINGEMENT AGAINST SOMATIK LABS AND ITS SUBSIDIARIES AND AFFILIATES OR THESOMATIK LABS SERVICES, RELATING TO OR ARISING OUT OF YOUR BREACH OF ANY PROVISION OF THESE TERMS, YOUR MISUSE OF THE SOMATIK LABS OR SOMATIK LABS SERVICES, OR YOUR UNAUTHORIZED MODIFICATION OR ALTERATION OF PRODUCTS OR SOFTWARE.

14.2 WARRANTY AND WARRANTY DISCLAIMER

SOMATIK LABS HAS A LIMITED WARRANTY ON WHEREBY SOMATIK LABS WARRANTS TO AND ONLY TO YOU THAT THE ON SOMATIK LABS MEDICAL DIAGNOSTIC AND UV-C SANITISATION DEVICES, INCLUDING TEMPERATURE THERMOMETERS, BLOOD PRESSURE MONITORS, UV-C SANITISERS AND UV-C SANITISING WANDS/CASES/TRAVEL BAGS WILL BE FREE FROM DEFECTS IN MATERIALS AND WORKMANSHIP FOR ONE (1) YEAR FROM THE DATE OF YOUR PURCHASE OF THE ON SOMATIK LABS MEDICAL DIAGNOSTIC AND UV-C SANITISATION DEVICES, INCLUDING TEMPERATURE THERMOMETERS, BLOOD PRESSURE MONITORS, UV-C SANITISERS AND UV-C SANITISING WANDS/CASES/TRAVEL BAGS (UNLESS A LONGER WARRANTY PERIOD IS REQUIRED BY LAW). THE SPECIFICS OF THIS SOMATIK LABS LIMITED WARRANTY ARE AT [HTTPS://WWW.SOMATIKLABS.COM/WARRANTY](https://www.somatiklabs.com/warranty). TO THE EXTENT POSSIBLE UNDER GOVERNING LAW, OTHER THAN THE ABOVE PRODUCT WARRANTY FOR THE SOMATIK LABS MEDICAL DIAGNOSTIC AND UV-C SANITISATION DEVICES,

INCLUDING TEMPERATURE THERMOMETERS, BLOOD PRESSURE MONITORS, UV-C SANITISERS AND UV-C SANITISING WANDS/CASES/TRAVEL BAGS, YOU UNDERSTAND AND AGREE THAT THE SOMATIK LABS SERVICES ARE PROVIDED ON AN "AS IS" AND "AS AVAILABLE" BASIS. SOMATIK LABS MAKES NO WARRANTY THAT THE SOMATIK LABS PRODUCT AND SERVICES WILL MEET YOUR REQUIREMENTS OR THAT USE OF THE SOMATIK LABS SERVICES WILL BE UNINTERRUPTED, TIMELY, SECURE OR ERROR-FREE NOR DOES SOMATIK LABS MAKE ANY WARRANTY AS TO THE ACCURACY OR RELIABILITY OF ANY INFORMATION OBTAINED THROUGH THE SOMATIK LABS SERVICES (INCLUDING THIRD PARTY CONTENT), THAT ANY DEFECTS IN THE SOMATIK LABS SERVICES WILL BE CORRECTED OR THAT THE SOMATIK LABS PRODUCTS OR SOMATIK LABS SERVICES WILL BE COMPATIBLE WITH ANY OTHER SPECIFIC HARDWARE OR SERVICE. FURTHER, SOMATIK LABS DOES NOT WARRANT THAT THE SOMATIK LABS SERVICES OR THE SOMATIK LABS SERVERS THAT PROVIDE YOU WITH DATA AND CONTENT ARE FREE OF VIRUSES OR OTHER HARMFUL COMPONENTS. SOMATIK LABS ALSO ASSUMES NO RESPONSIBILITY FOR AND SHALL NOT BE LIABLE FOR ANY DAMAGES CAUSED BY VIRUSES THAT MAY INFECT YOUR SOMATIK LABS MEDICAL DIAGNOSTIC AND UV-C SANITISATION DEVICES, INCLUDING TEMPERATURE THERMOMETERS, BLOOD PRESSURE MONITORS, UV-C SANITISERS AND UV-C SANITISING WANDS/CASES/TRAVEL BAGS, COMPUTER SOFTWARE, OR OTHER HARDWARE.


YOU AGREE THAT SOMATIK LABS IS NOT A PROVIDER OF MEDICAL SERVICES AND THAT SOMATIK LABS IS NOT PROVIDING ME WITH MEDICAL DIAGNOSIS OF ANY TYPE. IN THE EVENT OF ANY LOSS, DAMAGE OR INJURY, YOU WILL NOT LOOK TO SOMATIK LABS TO COMPENSATE YOU OR ANYONE ELSE. YOU RELEASE AND WAIVE FOR YOURSELF AND YOUR INSURER ALL

SUBROGATION AND OTHER RIGHTS TO RECOVER AGAINST SOMATIK LABS ARISING AS A RESULT OF THE PAYMENT OF ANY CLAIM FOR LOSS, DAMAGE OR INJURY.

SOMATIK LABS's EQUIPMENT AND SERVICES DO NOT CAUSE AND CANNOT ELIMINATE OCCURRENCES OF CERTAIN EVENTS, INCLUDING, BUT NOT LIMITED TO MEDICAL PROBLEMS OR MEDICAL INCIDENTS. SOMATIK LABS MAKES NO GUARANTEE OR WARRANTY, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, THAT THE SOMATIK LABS EQUIPMENT AND SERVICES PROVIDED WILL DETECT OR AVERT SUCH INCIDENTS OR THEIR CONSEQUENCES. SOMATIK LABS DOES NOT UNDERTAKE ANY RISK THAT YOU OR PROPERTY, OR THE PERSON OR PROPERTY OF OTHERS, MAY BE SUBJECT TO INJURY OR LOSS IF SUCH AN EVENT OCCURS. THE ALLOCATION OF SUCH RISK REMAINS WITH YOU, NOT SOMATIK LABS.

OTHER THAN THE ABOVE PRODUCT WARRANTY FOR THE SOMATIK LABS MEDICAL DIAGNOSTIC AND UV-C SANITISATION DEVICES, INCLUDING TEMPERATURE THERMOMETERS, BLOOD PRESSURE MONITORS, UV-C SANITISERS AND UV-C SANITISING WANDS/CASES/TRAVEL BAGS, SOMATIK LABS AND ITS SUPPLIERS DISCLAIM ALL WARRANTIES OF ANY KIND, WHETHER EXPRESS, IMPLIED, OR STATUTORY, REGARDING THE SOMATIK LABS SERVICES, INCLUDING ANY IMPLIED WARRANTY OF TITLE, MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NONINFRINGEMENT OF THIRD PARTY RIGHTS. BECAUSE SOME JURISDICTIONS DO NOT PERMIT THE EXCLUSION OF IMPLIED WARRANTIES, THE LAST SENTENCE OF THIS SECTION MAY NOT APPLY TO YOU. SOMATIK LABS HEREBY FURTHER EXPRESSLY DISCLAIMS ALL LIABILITY FOR ANY CLAIMS FOR SERVICE FAILURES THAT ARE DUE TO NORMAL PRODUCT WEAR, PRODUCT MISUSE, ABUSE, PRODUCT

MODIFICATION, IMPROPER PRODUCT SELECTION OR YOUR NON-COMPLIANCE WITH ANY AND ALL APPLICABLE FEDERAL, STATE OR LOCAL LAWS. THIS WARRANTY AND WARRANTY DISCLAIMER GIVE YOU SPECIFIC LEGAL RIGHTS, AND YOU MAY HAVE OTHER RIGHTS THAT VARY BY STATE, PROVINCE, OR COUNTRY. OTHER THAN AS PERMITTED BY LAW, SOMATIK LABS DOES NOT EXCLUDE, LIMIT OR SUSPEND OTHER RIGHTS YOU HAVE, INCLUDING THOSE THAT MAY ARISE FROM THE NONCONFORMITY OF A SALES CONTRACT. FOR A FULL UNDERSTANDING OF YOUR RIGHTS YOU SHOULD CONSULT THE LAWS OF YOUR STATE, PROVINCE, OR COUNTRY. FOR OUR AUSTRALIAN CUSTOMERS: PLEASE NOTE THAT THIS WARRANTY IS IN ADDITION TO ANY STATUTORY RIGHTS IN AUSTRALIA IN RELATION TO YOUR GOODS WHICH, PURSUANT TO THE AUSTRALIAN CONSUMER LAW, CANNOT BE EXCLUDED.



Engineered for your well-being.

DIRECTED.

All Images, Information and Content are Copyright
© 2020 Directed Electronics Australia Pty Ltd.
All Rights Reserved.

www.directed.com.au | www.directed.co.nz



www.somatiklabs.com